


## DRAWER DIVIDER KIT DKA-3003

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

### DKA-3003 Divider Kit for 30" wide “A” drawer

Configuration: 12 spaces (6.2"wide), 6 spaces (13"wide)



**Dividers:** (stainless steel)  
10 (6.677" x 1.75")  
5 (13.47" x 1.75")

**Rails:** (finished in same manner as drawer)  
1 Left Side  
1 Right Side  
2 Centers

**Drawer Mats:**  
2 (6.677" wide)  
1 (13.47" wide)