


DRAWER DIVIDER KIT DKA-3603

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

DKA-3603 Divider Kit for 36" wide “A” drawer

Configuration: 12 spaces (7.7" wide), 6 spaces (16" wide)



Dividers: (stainless steel)
10 (8.177" x 1.75")
5 (16.47" x 1.75")

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
2 (8.177" wide)
1 (16.47" wide)