

## BENCH HIGH BASE CABINET WL-3011

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.



### WL-3011 FEATURES

**Six "A" drawer, one "B" drawer base cabinet.**

Cabinet Dimensions: 30" wide, 35 1/2" high, 22" deep

#### "A" Drawer

Drawer Face Dimensions: 28 3/8" wide, 3 1/4" high, 1" thick.  
Drawer Body Inside Dimensions: 25 7/8" w, 2" h, 19 1/2" deep.  
Weight Capacity: 75 lbs.

#### "B" Drawer

Drawer Face Dimensions: 28 3/8" wide, 6 1/16" high, 1" thick.  
Drawer Body Inside Dimensions: 25 7/8" w, 5" h, 19 1/2" deep.  
Weight Capacity: 75 lbs.

