

ISLAND BASE CABINET WL-5822-1

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

WL-5822 FEATURES

Eight "B" drawer, one "D" door base cabinet.

Cabinet Dimensions: 58" wide, 35 1/2" high, 22" deep

NOTE:

Island Cabinets are also available with standard full-width bases for use when not positioned at the end of an island.

"B" Drawer

Drawer Face Dimensions: 16 3/8" wide, 6 11/16" high, 1" thick.
Drawer Body Inside Dims.: 13 7/8" wide, 5" high, 19 1/2" deep.
Weight Capacity: 75 lbs.

"D" Door

Door Face Dimensions: 16 3/8" wide, 27 3/16" high, 1" thick.

