


DRAWER DIVIDER KIT DKB-2203

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

DKB-2203 Divider Kit for 22” wide “B” drawer

Configuration: 30 spaces (5.8” wide)



Dividers: (stainless steel)
27 (6.275” x 3.5”)

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
3 (6.275” wide)