


DRAWER DIVIDER KIT WDKA-2404

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKA-2404 Divider Kit for 24” wide “A” drawer

Configuration: 12 spaces (6” wide)



Dividers: (Stainless Steel)
9 (6.441" x 1.75")

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
3 (6.441" wide)