



DRAWER DIVIDER KIT WDKA-4202

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKA-4202 Divider Kit for 42” wide “A” drawer

Configuration: 12 spaces (12.1” wide)



Dividers: (Stainless Steel)
9 (12.557" x 1.75")

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
3 (12.557"wide)