


DRAWER DIVIDER KIT WDKB-4802

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKB-4802 Divider Kit for 48” wide “B” drawer

Configuration: 12 spaces (14.1” wide)



Dividers: (Stainless Steel)
9 (14.557” x 3.5”)

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
3 (14.557” wide)

