



DRAWER DIVIDER KIT WDKB-4803

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKB-4803 Divider Kit for 48” wide “B” drawer

Configuration: 8 spaces (10.4” wide), 4 spaces (21.3” wide)



Dividers: (Stainless Steel)
6 (10.802” x 3.5”)
3 (21.72” x 3.125”)

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
2 (10.802” wide)
1 (21.72” wide)