


CHANNEL CABINET SPACERS

Channel Cabinet Spacers are available in two heights to correspond with the two available heights of standard Base Cabinets. Placing the Channel Cabinet Spacers between the Base Cabinets provides a flush filler panel. Each Channel Cabinet Spacer is finished in the same manner as the Base Cabinets.

Channel Cabinet Spacers for laboratory style (wetlab) cabinets are available in two heights:

WCCS-15 for bench-height cabinets


WCCS-10 desk-height cabinets



BENCH-HIGH (WCCS-15) Dimensions

variable width

3" top 2"




Bottom 6 1/2" painted black to match cabinet base

DESK-HIGH (WCCS-10) Dimensions

variable width

3" top 2"



Bottom 2 1/2" painted black to match cabinet base