

**MODESTY PANELS**

Each Modesty Panel is a reinforced assembly and is finished in the same manner as the Base Cabinets. Modesty Panels extend between Base Cabinets or Panel Legs to provide stability to some workbenches. Modesty Panels also provide some privacy and can be used to hide plumbing, wiring or duct work behind the workbench.


**Modesty Panels** for laboratory style (wetlab) cabinets are available in four heights:

**WMP-15** for bench-height cabinets


**WMPA-15** bench-height cabinets with an Apron Drawer.

**WMP-10** desk-height cabinets

**WMPA-10** desk-height cabinets with an Apron Drawer.



**BENCH-HIGH (WMP-15)  
Modesty Panel Dimensions**



**DESK-HIGH (WMP-10)  
Modesty Panel Dimensions**

