

ACID STORAGE BASE CABINET - (ASC)

Acid Storage Cabinets safely store and vent acids and other corrosive liquids.

FEATURES:

- Durable epoxy-coated steel construction with corrosion-resistant polyethylene-lined interior
- Attractive glacier white exterior complements laboratory casework.
- Manual-closing, non-locking door(s)
- Include vacuum-formed PVC liner tray.
- Include four leveling feet and one 8" filler panel to increase cabinet depth from 22" to 30".
- Support loads up to 800 pounds.

OPTIONS:

- Vent Kit
- Shelf Kits
- Filler Panels
- ADA Height

REQUIREMENTS:

- Work Surface



SPECIFICATIONS:

18" WIDTH (LEFT-HINGE OR RIGHT-HINGE SINGLE DOOR)

ATTRIBUTES:

Cabinet Width: 18" (46 cm)

Door: Single, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Acid storage

SPECIFICATIONS:

Estimated Shipping Weight: 100.0 lbs

Estimated Shipping Weight metric: 45.4 kg

Dimensions: 18.0" w x 22.0" d

Dimensions metric: 45.7 x 55.9 cm

24" WIDTH (LEFT-HINGE OR RIGHT-HINGE SINGLE DOOR)

ATTRIBUTES:

Cabinet Width: 24" (61 cm)

Door: Single, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Acid storage

SPECIFICATIONS:

Estimated Shipping Weight: 110.0 lbs

Estimated Shipping Weight metric: 50.0 kg

Dimensions: 24.0" w x 22.0" d

Dimensions metric: 61.0 x 55.9 cm

30" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 30" (76 cm)

Door: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Acid storage

SPECIFICATIONS:

Estimated Shipping Weight: 145.0 lbs

Estimated Shipping Weight metric: 65.8 kg

Dimensions: 30.0" w x 22.0" d

Dimensions metric: 76.2 x 55.9 cm

36" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 36" (76 cm)

Door: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Acid storage

SPECIFICATIONS:

Estimated Shipping Weight: 155.0 lbs

Estimated Shipping Weight metric: 70.3 kg

Dimensions: 36.0" w x 22.0" d

Dimensions metric: 91.4 x 55.9 cm

48" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 48" (122 cm)

Door: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Acid storage

SPECIFICATIONS:

Estimated Shipping Weight: 170.0 lbs

Estimated Shipping Weight metric: 77.1 kg

Dimensions: 48.0" w x 22.0" d

Dimensions metric: 121.1 x 55.9 cm