

SOLVENT STORAGE BASE CABINET - (FSC)

Solvent Storage Cabinets safely store and vent solvents and other flammable liquids. In case of a building fire, these cabinets protect the contents inside them.

STANDARD FEATURES:

- Durable dual wall epoxy-coated steel construction with 1.5" air gap
- Attractive glacier white exterior complements laboratory casework.
- Include epoxy-coated steel shelf
- Include two rear vent connections with flame arrestors and closure plugs
- Include four leveling feet and one 8" filler panel to increase cabinet depth from 22" to 30".
- Support loads up to 800 pounds.
- Manual-closing door(s) with 3-position locking handle

OPTIONAL FEATURES :

- Self-closing door(s) with 3-position locking handle. Self-closing door(s) close automatically in case of fire when temperature reaches 74° C (165° F)
- Shelf Kits
- Filler Panels
- ADA Height

REQUIREMENTS:

- Work Surface



SPECIFICATIONS:

24" WIDTH (LEFT-HINGE OR RIGHT-HINGE SINGLE DOOR)

ATTRIBUTES:

Cabinet Width: 24" (61 cm)

Door: Single, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Solvent storage

SPECIFICATIONS:

Estimated Shipping Weight: 155.0 lbs

Estimated Shipping Weight metric: 59.0 kg

Dimensions: 24.0" w x 22.0" d

Dimensions metric: 61.0 x 55.9 cm

30" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 30" (76 cm)

Doors: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Solvent storage

SPECIFICATIONS:

Estimated Shipping Weight: 195.0 lbs

Estimated Shipping Weight metric: 70.0 kg

Dimensions: 30.0" w x 22.0" d

Dimensions metric: 76.2 x 55.9 cm

36" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 36" (91 cm)

Doors: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Solvent storage

SPECIFICATIONS:

Estimated Shipping Weight: 205.0 lbs

Estimated Shipping Weight metric: 82.0 kg

Dimensions: 36.0" w x 22.0" d

Dimensions metric: 91.4 x 55.9 cm

48" WIDTH (DOUBLE DOORS)

ATTRIBUTES:

Cabinet Width: 48" (122 cm)

Doors: Dual, Manual Close

Filler Panel Depth: 8.0" (20.3 cm)

Maximum Base Height: 36.75" (93.3 cm)

Minimum Base Height: 35.5" (90.2 cm)

Storage Type: Solvent storage

SPECIFICATIONS:

Estimated Shipping Weight: 250.0 lbs

Estimated Shipping Weight metric: 118.0 kg

Dimensions: 48.0" w x 22.0" d

Dimensions metric: 121.9 x 55.9 cm